

MOMentum

By Visual MED



Lets you deliver . . .
'Irresistible content'TM

. . . for FREE. Really.

Help your patients share this exciting time in their lives.
Send their ultrasound images, right to their smartphones.
MOMentum is Free for OBGYN practices.



Enhance your patient's practice experience.

Since the invention of ultrasound, expectant moms have been getting the first pictures of their new baby on little black and white thermal prints. The thermal print is 40 year old technology that starts to lose its quality as soon as it is removed from the printer. The thermal print is hard to scan and is hard to share with the important people in the lives of your patients.

MOMentum allows the OBGYN practice to send the 'all important first images' of your expectant mom's baby, digitally and securely, to their mobile device or computer. The app allows you to help your expectant moms **share this wonderful and important time** in their life with their friends and family, from the **convenience** of their smart phone. Their ultrasound images are anonymized and digital and they arrive in **full quality**.



The MOMentum app **eliminates** your need for thermal prints and a thermal printer.

You can even include your practice logo with the images. **Great free marketing.**

With the app, you can **text** or **email** your expectant moms their images. Enhancing their practice experience.

Once they have their images on their phones they can **download** their images in **full resolution and in perfect quality**.

From their smartphone's apps they can annotate their images and instantly share their images with the important people in their life using their favorite social media.

90 MOMentum video
https://youtu.be/KLUTAwgK_64

Set up is easy. All you need for set up, you probably already have, so you do not need to purchase any hardware:

- . DICOM ultrasound system
- . Network connection for your ultrasound system. LAN or direct connect.
- . PC or laptop to run the MOMentum app. WIN OS 7 / 10