

# Give yourself the EDGE in Imaging

The **EDGE 3 in 1** ultrasound system ..

- . Smaller than your iPhone.
- . Designed by providers and sonographers
- . Gives your diagnostic visualization on every patient encounter.
- . Boots up in 1 second.
- . Improves your patient confidence and compliance.
- . Enhances the value of your time and respects your patient's time.
- . Fits into your practice operating budget without expanding it.
- . Multi frequency Convex, Linear and Phased Array imaging formats
- . iOS and Droid operating system.
- . Download and updates are free at your App Store and Play Store.



Display	Phone or tablet screen
OS	iOS / Droid .. Free download and updates
MicroVue App	Operates as a standard phone app. Swipe left and right to change scan modes, gain, slide up and down to change depth ..
Probe Type	Convex array / Linear array / Sector
Scan Angle	Convex 60mm
Scan Length	Linear 40mm
Frequency	Convex probe 3.5 / 5.0 MHz, Linear probe 7.5 / 10MHz

Scan mode	B, B/M, Color, PW
Measure	Length, Area, Velocity, HR, S/D, OB
Puncture Assist	In-plane, out-of-plane, needle enhance
Frame Rate	20 fps
Cine play	10 Seconds
Power	Built-in lithium battery, 2h, work 2.5h
Charging	Wireless charger
Display Depth	Convex 90 ~ 305mm, Linear 20 ~ 100mm
Image Adjust	Gain, Depth, DR, Focus, SRI, Harmonic

Less than \$260.00 per month\*

Includes: MOMENTUM modality Worklist, Slim line mobile stand, on-line training and support and two-year warranty.

\*Credit approval.



Call for more information or a demo .. (800) 821-3230 [www.visual-med.com](http://www.visual-med.com) Click on EDGE 3 in 1